


55.GA UWAVS Optical Alignment

Ad Verlaan
Systems Engineer
TNO

Presented at the
UWAVS Preliminary Design Review
September 17-19, 2024



Objectives & requirement

UWAVS alignment


- **Optical alignment is key**
 - To achieve the imaging resolution
 - After initial alignment
 - And to accommodate vessel deformation
- **Optical alignment is performed**
 - Using alignment laser and targets
 - Using active alignment system
- **Optical alignment**
 - Is part of the overall alignment plan, feeds into assembly and qualification
 - Is considered in tolerance analysis
 - Accuracies are proposed

Key Assumptions

- **Alignment**
 - DSM interfaces errors well known & small
 - Laser tracker available in PIF and port cell
- **Active alignment**
 - In-vessel contrast in visible and IR
 - Vessel deformation shows smooth behavior
 - Vibrations are not corrected

IDM #	Reference Documents
AQLWZW	55.GA UWAVS-35 Optical Model Results
AK7KBE	55.GA UWAVS-35 Onsite Assembly Plan
WL7A5Y	55.GA UWAVS-35 Onsite Testing and Commissioning Plan

System initial alignment and role of optical alignment



IDM #	Reference Documents
AQLWZW	55.GA UWAVS-35 Optical Model Results
AK7KBE	55.GA UWAVS-35 Onsite Assembly Plan
WL7A5Y	55.GA UWAVS-35 Onsite Testing and Commissioning Plan

- **Mechanical alignments**
 - Measure & correct shims or adjust element
 - Mainly elements & modules
 - CMM measurements
- **Laser tracker alignment**
 - Modules in EPP, ISS & PCSS
 - Mainly larger units
- **Optical alignments**
 - A40 pre-alignment
 - A50 mirror alignment
 - A60 pre-alignment
 - A80, A90 & A100 alignment
 - To pre-align compensator and
 - Achieve better element alignment
- **Optical verification**
 - Optical axes
 - Imaging performance

A50 Alignment to improve IOT imaging

IDM #	Reference Documents
AQLWZW	55.GA UWAVS-35 Optical Model Results
AK7KBE	55.GA UWAVS-35 Onsite Assembly Plan
WL7A5Y	55.GA UWAVS-35 Onsite Testing and Commissioning Plan

- **Goals**


- Adjust A50 mirrors in 2 degrees of freedom (tip and tilt)
- To place optical axis centred on the mirrors
- To improve the imaging IOT imaging resolution

- **Approach**

- Set up alignment laser at nominal position & direction using laser tracker and fiducial
- Measure transmitted beam position & direction
- Iteratively adjust A40 tip-tilt & piston

- **Tools**

- Alignment laser with small beam
- Laser tracker to measure positions in 3D
- Fiducial to link beam to position
- Beam position markers for each mirror to mark the spot



Example optical alignment of A50 mirror group mirrors

A40 pre-alignment for ISS integration

IDM #	Reference Documents
AQLWZW	55.GA UWAVS-35 Optical Model Results
AK7KBE	55.GA UWAVS-35 Onsite Assembly Plan
WL7A5Y	55.GA UWAVS-35 Onsite Testing and Commissioning Plan

- **Goals**


- Adjust A40 mirror in 5 degrees of freedom
- To place optical axis right wrt mechanical axis
- In preparation to ISS integration and verification

- **Approach**

- Set up alignment laser at nominal position & direction using laser tracker and fiducial
- Measure transmitted beam position & direction
- Iteratively adjust A40 tip-tilt & piston

- **Tools**

- Alignment laser with small beam
- Laser tracker to measure positions in 3D
- Fiducial to link beam to position



A60 pre-alignment for ISS integration

IDM #	Reference Documents
AQLWZW	55.GA UWAVS-35 Optical Model Results
AK7KBE	55.GA UWAVS-35 Onsite Assembly Plan
WL7A5Y	55.GA UWAVS-35 Onsite Testing and Commissioning Plan

- **Goals**


- Adjust A60 mirror in 2 times 2 degrees of freedom (tip-tilt)
- To place optical axis out wrt the axis in
- In preparation to ISS integration and verification

- **Approach**

- Set up alignment laser at nominal position & direction using laser tracker and fiducial
- Measure transmitted beam position & direction
- Iteratively adjust A40 tip-tilt & piston

- **Tools**

- Alignment laser with small beam
- Laser tracker to measure positions in 3D
- Fiducial to link beam to position



BEOC mirror & IR breadboard alignment for integration

- **Goals**

- Adjust A80, A100 mirrors and IR baseplate in 2 degrees of freedom
- To place optical axis right on the imagers
- In preparation to integration & verification



- **Approach**

- Set up alignment laser at nominal position & direction using laser tracker and fiducial
- Measure transmitted beam position & direction
- Iteratively adjust A40 tip-tilt & piston

- **Tools**


- Alignment laser with small beam
- Fiducials and reference tool to mark beam position

IDM #	Reference Documents
AQLWZW	55.GA UWAVS-35 Optical Model Results
AK7KBE	55.GA UWAVS-35 Onsite Assembly Plan
WL7A5Y	55.GA UWAVS-35 Onsite Testing and Commissioning Plan



A90 and A110 lens centration to improve imaging

- **Goal**
 - Reduce positioning tolerances of lenses in the lens barrels
 - To achieve imaging resolution
- **Approach**
 - Align barrel
 - Align lens
 - Lock position & move to next lens
- **Tool**
 - Lens centration bench



IDM #	Reference Documents
AQLWZW	55.GA UWAVS-35 Optical Model Results
AK7KBE	55.GA UWAVS-35 Onsite Assembly Plan
WL7A5Y	55.GA UWAVS-35 Onsite Testing and Commissioning Plan

Active alignment system overview

- **System needs**

- Correct beam position, pointing and focus
- For initial alignment
- And to compensate for vessel deformation



- **Implementation**

- A40-M1: Tip-tilt
- A40-M2: Tip-tilt & piston
- Visible camera focus adjustment

- **Operation**

- Bandwidth $\sim 1\text{Hz}$, so only thermal deformations
- Use pre-determined thermal correction
- And image position on detectors

Alignment system implementation




Alignment system Operation

IDM #	Reference Documents
AQLWZW	55.GA UWAVS-35 Optical Model Results

Performance evaluation: A40 actuation range

- **Actuation range is determined by adding:**
 - Tolerance analysis alignment ranges
 - STOP deformation correction, this is in one direction only
- **Outcome**
 - Tip-tilt ranges are fine
 - Focus requires non-central use of range (shift center of range)



Overview	A40-M1		A40-M2		
Coordinates	ZLCS	ZLCS	TGCS Port 0	ZLCS	ZLCS
	Rx [mrad]	Ry [mrad]	dz [mm]	Rx [mrad]	Ry [mrad]
STOP IOT-Zemax temp analysis (52°C) BEOC-Zemax temp analysis (37°C) Vacuum vessel deformation	-0.31	-23.35	-19.20	0.23	-21.80
Tolerances (-) direction	-25.20	-32.37	-4.70	-31.21	-38.57
Tolerances (+) direction	25.20	32.37	4.70	31.21	38.57
Total range required: Thermal + tolerances	50.71	88.10	28.60	62.65	98.94

Active alignment system has sufficient range

IDM #	Reference Documents
AQLWZW	55.GA UWAVS-35 Optical Model Results

Performance evaluation: A40 Actuator resolution

- **Actuation resolution need**
 - Sub-pixel alignment on detector
 - To minimize impact on data interpretation
- **Implemented resolution**
 - $\pm 1\mu\text{m}$ on focus
 - $\pm 0.65\text{mrad}$ on tip-tilt
- **Outcome**
 - Better angular resolution required, foreseen for FD-phase

Active alignment resolution needs improvement

IDM #	Reference Documents
AQLWZW	55.GA UWAVS-35 Optical Model Results




Image displacement by A40 mirror tilt

Targeted resolutions for A40 actuators	
Tilt actuators	<0,008 mrad
Piston	<0,045 mm

Targeted actuator resolution


Performance evaluation: A40 Mirror size

- **Beam shift is determined by adding**
 - Thermal deformation beam shift (From STOP analysis)
 - Initial alignment beam shift (tolerance analysis)
 - Pivot point effect (rotation point is 69.5mm behind surface)
- **Outcome**
 - Smallest margin is ~5mm
 - No vignetting will occur on A40 mirrors during alignment and operation


Active alignment mirrors sufficiently large

IDM #	Reference Documents
AQLWZW	55.GA UWAVS-35 Optical Model Results


Nominal footprint on A40-M1




Max. footprint on A40-M1 after tolerances



Nominal footprint on A40-M2



Max. footprint on A40-M2 after tolerances



Camera focus adjustment visible cameras

- **Focus correction difference visible and IR**
 - A40 adjusts all channels at same time, yet
 - Thermal behaviour Visible and IR channels BEOC different
- **Implementation**
 - Also visible camera focus added
 - Use A40 for general focus
 - And visible camera for difference only
 - Use model based LUT and BEOC temperature for adjustment

Active alignment camera focusing implemented

IDM #	Reference Documents
AQLWZW	55.GA UWAVS-35 Optical Model Results

Alignment	Vis-1	IR-1
A40-M2 focus	2.4	-1.06
A40-M1 X tilt	0.285	-0.15
A40-M1 Y tilt	0.003	0.02
A40-M2 X tilt	0.285	-0.14
A40-M2 Y tilt	-0.003	0.03
Camera focus	-0.352	0

Final system compensator range for IR-1 and Vis-1 channels

Active alignment system operation

- **Active alignment system**

- Combines thermal and initial alignment
- During operation only has image position data
- Focus correction requires some scanning
- Angular correction


- **Approach**

- Perform initial alignment scan & optimization of all degrees of freedom after installation (T=20DegC)
- Repeat alignment scan & optimization of all degrees of freedom on warm vessel (T=70DegC)
- Create LUT for temperature correction
- Add image position correction during operation/between pulses

- **Open issues**

- Available contrast in-vessel around pulses

Active alignment process flow



Active alignment process flow defined

IDM #	Reference Documents
AQLWZW	55.GA UWAVS-35 Optical Model Results

Summary slide

UWAVS alignment

- **System optical alignment is presented for**
 - Initial system alignment
 - Active alignment system for vessel deformation
 - Using two flat mirror tip, tilt and piston adjustment
- **Initial system optical alignment**
 - Helps to achieve accurate optics alignment
 - Prepares for system qualification measurements
- **Active alignment system**
 - Has sufficient actuator range
 - Mirrors are sufficiently large
 - Actuator resolution will be improved in FD-phase

Open Issues & way forward

- **Alignment**
 - DSM interfaces errors well known & small
 - Access for Laser tracker measurement in PIF and port cell
 - Plan: discuss with PI in FD-phase
- **Active alignment**
 - Actuator resolution to be improved (in FD-phase)
 - Plan: Update design in FD-phase
 - Illumination & contrast around pulses for alignment
 - Plan: Inventory options in FD-phase